



















**3**

Assay Data Sheet - Expected Ranges  
 BioProfile® FLEX Controls - LEVEL 3

**3**

Lot Number:  610351		Exp. Date: 2017-10	
		Expected Range	
Parameter	Units	Lower Limit	Upper Limit
pH	-	 7.306	 7.456
PO <sub>2</sub>	mm Hg	 132.3	 175.1
PCO <sub>2</sub>	mm Hg	 12.0	 21.6
Gluc	g/L	 11.62	 18.38
Lac	g/L	 3.87	 6.13
NH <sub>4</sub> <sup>+</sup>	mmol/L	 7.49	 11.99
Na <sup>+</sup>	mmol/L		
K <sup>+</sup>	mmol/L		
Ca <sup>++</sup>	mmol/L		
Osmolality	mOsm/kg H <sub>2</sub> O	 605	 621

*nova*  
 biomedical

200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800

LFN 08241A 2007-08