


































**1**

## Assay Data Sheet - Expected Ranges BioProfile® FLEX Controls - LEVEL 1

**1**

Lot Number:  105621				Exp. Date: 2012-11	
Expected Range				+/- 3 SD Expected Range	
Parameter	Units	Lower Limit	Upper Limit	Lower Limit	Upper Limit
pH	-	 6.615	 6.715	 6.590	 6.740
PO <sub>2</sub>	mm Hg	 59.6	 72.0	 56.5	 75.1
PCO <sub>2</sub>	mm Hg	 59.8	 73.4	 56.4	 76.8
Gluc	g/L	 0.75	 1.25	 0.62	 1.38
Lac	g/L	 0.78	 1.28	 0.65	 1.41
NH <sub>4</sub> <sup>+</sup>	mmol/L	 0.61	 1.11	 0.48	 1.24
Na <sup>+</sup>	mmol/L	 69.7	 85.3		
K <sup>+</sup>	mmol/L	 1.71	 2.29		
Ca <sup>++</sup>	mmol/L	 0.50	 0.74		
Osmolality	mOsm/kg H <sub>2</sub> O	 214	 224		

*nova*  
biomedical

200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800

LFN 08239A 2007-08