

1

Assay Data Sheet - Expected Ranges  
BioProfile® FLEX Controls - LEVEL 1

1

Lot Number:



Exp. Date: 2021-01

Parameter	Units	Expected Range		+/- 3 SD Expected Range	
		Lower Limit	Upper Limit	Lower Limit	Upper Limit
pH	-	 6.626	 6.726	 6.601	 6.751
PO <sub>2</sub>	mm Hg	 56.5	 68.9	 53.4	 72.0
PCO <sub>2</sub>	mm Hg	 58.2	 71.8	 54.8	 75.2
Gluc	g/L	 0.75	 1.25	 0.62	 1.38
Lac	g/L	 0.78	 1.28	 0.65	 1.41
NH <sub>4</sub> <sup>+</sup>	mmol/L	 0.61	 1.11	 0.48	 1.24
Na <sup>+</sup>	mmol/L	 69.7	 85.3		
K <sup>+</sup>	mmol/L	 1.71	 2.29		
Ca <sup>++</sup>	mmol/L	 0.50	 0.74		
PO <sub>4</sub>	mmol/L	 3.42	 4.62		
p	mg/L	 106.0	 143.2		
Osmolality	mOsm/kg H <sub>2</sub> O	 217	 227		