

5

Assay Data Sheet - Expected Ranges  
BioProfile® FLEX Controls - LEVEL 5

5

Lot Number:  832437		Exp. Date: 2020-05	
Parameter	Units	Expected Range	
		Lower Limit	Upper Limit
Gln	mmol/L	 4.39	 5.89
Glu	mmol/L	 1.80	 2.40
Gluc	g/L	 8.11	 12.11
Lac	g/L	 4.46	 5.96



200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800