


























2

Assay Data Sheet - Expected Ranges
BioProfile® FLEX Controls - LEVEL 2

2

Lot Number:  509015		Exp. Date: 2016-09	
Parameter	Units	Expected Range	
		Lower Limit	Upper Limit
pH	-	 6.890	 6.990
PO ₂	mm Hg	 102.5	 122.5
PCO ₂	mm Hg	 41.4	 50.2
Gluc	g/L	 4.25	 5.75
Lac	g/L	 2.55	 3.45
NH ₄ ⁺	mmol/L	 3.85	 5.65
Na ⁺	mmol/L	 122.5	 140.9
K ⁺	mmol/L	 3.58	 4.42
Ca ⁺⁺	mmol/L	 0.67	 0.93
PO ₄	mmol/L	 0.75	 1.15
P	mg/L	 23.3	 35.7
Osmolality	mOsm/kg H ₂ O	 395	 407

Lower Limit

+/- 3 SD Expected Range

Upper Limit

LFN 08240B 2011-03

200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800

