

1

Assay Data Sheet - Expected Ranges
BioProfile® FLEX Controls - LEVEL 1

1

Lot Number:



Exp. Date: 2023-02

Parameter	Units	Expected Range		+/- 3 SD Expected Range	
		Lower Limit	Upper Limit	Lower Limit	Upper Limit
pH	-	 6.628	 6.728	 6.603	 6.753
PO ₂	mm Hg	 54.6	 67.0	 51.5	 70.1
PCO ₂	mm Hg	 59.0	 72.6	 55.6	 76.0
Gluc	g/L	 0.75	 1.25	 0.62	 1.38
Lac	g/L	 0.78	 1.28	 0.65	 1.41
NH ₄ ⁺	mmol/L	 0.61	 1.11	 0.48	 1.24
Na ⁺	mmol/L	 69.7	 85.3		
K ⁺	mmol/L	 1.71	 2.29		
Ca ⁺⁺	mmol/L	 0.50	 0.74		
PO ₄	mmol/L	 3.44	 4.64		
p	mg/L	 106.7	 143.9		
Osmolality	mOsm/kg H ₂ O	 217	 227		