





































1

Assay Data Sheet - Expected Ranges
BioProfile® FLEX Controls - LEVEL 1

1

Lot Number: 
703767

Exp. Date: 2018-08

Parameter	Units	Expected Range		+/- 3 SD Expected Range	
		Lower Limit	Upper Limit	Lower Limit	Upper Limit
pH	-	 6.651	 6.751	 6.626	 6.776
PO ₂	mm Hg	 58.8	 71.2	 55.7	 74.3
PCO ₂	mm Hg	 55.7	 69.3	 52.3	 72.7
Gluc	g/L	 0.75	 1.25	 0.62	 1.38
Lac	g/L	 0.78	 1.28	 0.65	 1.41
NH ₄ ⁺	mmol/L	 0.61	 1.11	 0.48	 1.24
Na ⁺	mmol/L	 69.7	 85.3		
K ⁺	mmol/L	 1.71	 2.29		
Ca ⁺⁺	mmol/L	 0.50	 0.74		
PO ₄	mmol/L	 3.62	 4.82		
p	mg/L	 112.2	 149.4		
Osmolality	mOsm/kg H ₂ O	 218	 228		



200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800