






































2

**Assay Data Sheet - Expected Ranges
BioProfile® FLEX Controls - LEVEL 2**

2

Lot Number:  706637		Exp. Date: 2018-09	
Parameter	Units	Expected Range	
		Lower Limit	Upper Limit
pH	-	 6.907	 7.007
PO ₂	mm Hg	 98.6	 118.6
PCO ₂	mm Hg	 35.6	 44.4
Gluc	g/L	 4.25	 5.75
Lac	g/L	 2.55	 3.45
NH ₄ ⁺	mmol/L	 3.85	 5.65
Na ⁺	mmol/L	 122.5	 140.9
K ⁺	mmol/L	 3.58	 4.42
Ca ⁺⁺	mmol/L	 0.67	 0.93
PO ₄	mmol/L	 0.76	 1.16
P	mg/L	 23.6	 36.0
Osmolality	mOsm/kg H ₂ O	 396	 408
		+/- 3 SD	Expected Range
		Lower Limit	Upper Limit
		 6.882	 7.032
		 93.6	 123.6
		 33.5	 46.5
		 3.87	 6.13
		 2.32	 3.68
		 3.35	 6.15



200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800

LFN 08240B 2011-03