


































3

Assay Data Sheet - Expected Ranges
BioProfile® FLEX Controls - LEVEL 3

3

Lot Number:  805868		Exp. Date: 2019-08	
Parameter	Units	Expected Range	
		Lower Limit	Upper Limit
pH	-	 7.309	 7.409
PO ₂	mm Hg	 133.3	 161.9
PCO ₂	mm Hg	 15.6	 22.0
Gluc	g/L	 12.75	 17.25
Lac	g/L	 4.25	 5.75
NH ₄ ⁺	mmol/L	 8.24	 11.24
Na ⁺	mmol/L	 166.7	 203.3
K ⁺	mmol/L	 5.55	 6.85
Ca ⁺⁺	mmol/L	 2.25	 2.85
Osmolality	mOsm/kg H ₂ O	 604	 620
		+/- 3 SD Expected Range	
		Lower Limit	Upper Limit
		 7.284	 7.434
		 126.2	 169.0
		 14.0	 23.6
		 11.62	 18.38
		 3.87	 6.13
		 7.49	 11.99



200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800