










5

Assay Data Sheet - Expected Ranges
BioProfile® FLEX Controls - LEVEL 5

5

Lot Number:  626740		Exp. Date: 2018-03	
Parameter	Units	Expected Range	
		Lower Limit	Upper Limit
Gln	mmol/L	 4.55	 6.05
Glu	mmol/L	 1.80	 2.40
Gluc	g/L	 8.29	 12.29
Lac	g/L	 4.40	 5.90

nova
biomedical

200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800