









5

Assay Data Sheet - Expected Ranges
BioProfile® FLEX Controls - LEVEL 5

5

Lot Number:  701057		Exp. Date: 2018-07	
Parameter	Units	Expected Range	
		Lower Limit	Upper Limit
Gln	mmol/L	 4.53	 6.03
Glu	mmol/L	 1.80	 2.40
Gluc	g/L	 8.19	 12.19
Lac	g/L	 4.44	 5.94

nova
biomedical

200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800

LFN 08243B 2008-01